




































**Medical Temperature Probe**

Part No	Description	Reference Picture	Compatibility	Original part number
PTP-10001A	YSI400, 25°C @2.2 52KΩ, single thermistor ;	 adult skin surface probe	Philips A3 series patient monitor; M3926A, M3927A, M3928A, M3929A; Marquette 7000/7010, GE Eagle, Solar, Tram, Transcope, MAC-Lab Systems; GE part number: M1024254	409B
PTP-10001B	Reusable temperature probes; 6.3mm plug diameter; phone plug; GE part number: 21729402	 adult rectal/esophageal		401
PTP-10001C	Probe, rectal/esophageal, 4 mm diameter, 400 Series;	 pediatric skin surface probe		427
PTP-10001D	21730002 Probe, rectal/esophageal, 3 mm	 pediatric rectal/esophageal		402
PTP-10001E		 adapter cable, using with disposable temperature probe; cable material is TPU		








PTP-10002A	<b>YSI 700 Series Reusable temperature probes with Dual Thermistors; 25°C@6K/30KΩ</b>	 adult skin surface probe	Spacelabs: 90496 Module, compatible with Spacelabs P/N: 700-0031-00 Marquette 7000/7010, GE Eagle, Solar, Tram, Tramscope, MAC-Lab Systems, compatible with GEMS P/N:9445-001	709B
PTP-10002B		 adult rectal/esophageal		701
PTP-10002C		 pediatric skin surface probe		
PTP-10002D		 pediatric rectal/esophageal		
PTP-10003A	<b>10K Series Reusable temperature probes with Single Thermistors</b>	 adult skin surface probe	Suitable for China OEM patient monitor: Tianrong: TR900D/E; Anke: ASC553A3, ASC553; Comen: star 8000A/B/C, star 5000, 5000B/C	
PTP-10003B		 adult rectal/esophageal		
PTP-10003C		 pediatric skin surface probe		
PTP-10003D		 pediatric rectal/esophageal		

PTP-10004A		 <p>adult skin surface probe</p>	Creative patient monitor: PC-1000, 5000, 9000, UP-6000,	
PTP-10004B	25°C@10 KΩ, Redel 2p	 <p>adult rectal/esophageal</p>		
PTP-10005A		adult skin surface probe		
PTP-10005B	Philips Reusable temperature probes with Single Thermistors;	adult rectal/esophageal		
PTP-10005C	21078A; 21075A; 21076A	 <p>pediatric skin surface probe</p>	<p>Temperature Module: M1029A, C3: 862474, 862478, SureSigns: 863064, 863065, 863066, 863068, 863086, 863087, 863088; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A</p>	
PTP-10005D		 <p>pediatric rectal/esophageal</p>		
PTP-10005E				








PTP-10005F			adapter cable using with disposable temperature probe, cable material is TPU;	
PTP-10006A	S&W Artema Reusable temperature probes, Single Thermistor	 <p>adult skin surface probe</p>	Diascope G2	
PTP-10006B		 <p>adult rectal/esophageal</p>		
PTP-10006C		 <p>pediatric skin surface</p>		
PTP-10006D		 <p>pediatric rectal/esophageal</p>		
PTP-10006E			adapter cable, using with disposable temperature probe; cable material is TPU	
PTP-10007A		<p>adult skin surface probe</p> 		7498939






PTP-10007B		 <p>www.panswmed.com</p>		
PTP-10007C	25°C@2.25KΩ	 <p>adult rectal/esophageal</p>	Siemens medical SC 6002XL, SC7000, SC8000, SC9000XL	
PTP-10007D		 <p>pediatric skin surface</p>		
PTP-10007E		 <p>pediatric rectal/esophageal</p>	adapter cable, using with disposable temperature probe; cable material is TPU	
PTP-10007F	5198333		adapter cable, to use YSI400 temperature probe on Drager machines	
PTP-10008A	Air Shields 25°C@6.867KΩ	 <p>adult skin surface probe</p>	25°C@6.867KΩ, Din3p C-100	
PTP-10008B		 <p>adult rectal/esophageal</p>		





PTP-10009A	15K Series Reusable temperature probes with	 <p>Pan-SW adult skin surface probe</p>		
PTP-10009B	Single Thermistors	 <p>Pan-SW adult rectal/esophageal</p>		
PTP-10009C	50K Series Reusable temperature probes with	 <p>Pan-SW adult skin surface probe</p>		
PTP-10009D	Single Thermistors	 <p>Pan-SW adult rectal/esophageal</p>		
PTP-10010A	Mindray PM6800 series patient monitors	 <p>Pan-SW adult skin surface probe</p>	PM6800 Series For BeneView, iPM, iMEC series monitors, BeneHeart defibrillator	0011-30-37393 MR403B
PTP-10010B		 <p>Pan-SW adult rectal/esophageal</p>		0011-30-37392 MR401B
PTP-10010C		 <p>Pan-SW pediatric skin surface</p>		0011-30-37395 MR404B











PTP-10010D		 <p>Pan-SW</p>		0011-30-37394 MR402B
PTP-10010E		 <p>Pan-SW</p>		adapter cable, using with disposable temperature probe; cable material is TPU
PTP-10001A		 <p>Pan-SW</p>		0011-30-90442 MR403
PTP-10001B	Mindray, <b>YSI400</b> , 25°C @2.2 52KΩ, single thermistor	 <p>Pan-SW</p>	Mindray all module, except PM6800 Series; PM8000; PM9000; MEC series, PM series	0011-30-90440 MR401
PTP-10001C	; Reusable temperature probes	 <p>Pan-SW</p>		0011-30-90443 MR404
PTP-10001D		 <p>Pan-SW</p>		0011-30-90441 MR402
PTP-10001E		 <p>Pan-SW</p>		


PTP-10011A	GE YSI 400 compatible, integrated probe from machine to patient;	adult skin surface probe	Dash, Eagle, Solar and TRAM-Monitors; Datex Ohmeda S/5 E-modules;	
PTP-10011B		adult rectal/esophageal		
PTP-10011C		pediatric skin surface		
PTP-10011D		pediatric rectal/esophageal		
PTP-10011E		adapter cable, using with disposable temperature probe; cable material is TPU		
PTP-10011F		adapter cable, dual temperature probe		
PTP-10011G		adapter to YSI 400 temperature probe;		





PTP-10012A		adult skin surface probe		
PTP-10012B	GE YSI 700 series compatible; integrated probe from machine to patient;	adult rectal/esophageal	Dash, Eagle, Solar and TRAM-Monitors; Datex Ohmeda S/5 E-modules;	
PTP-10012C		pediatric skin surface		
PTP-10012D		pediatric rectal/esophageal		
PTP-10013A	Ohmeda incubator temperature probe; Plug diameter: 3.5mm;	<p style="text-align: center;">skin surface probe</p> 	<p>Compatible with: Incubators: Ohmeda/Ohio IC, GC Care Plus.Warmers: Ohmeda Panda Baby Warmers and Ohmeda 3000, 4000 &amp; 5000 Series of Infant Warmer Systems</p>	<p>208-0697-700, 6600-0628-700; LA003</p>
PTP-10014A		adult skin surface probe		

PTP-10014B	spacelabs new type temperatu re probe	adult rectal/esophageal	
PTP-10014C		pediatric skin surface	
PTP-10014D		pediatric rectal/esophageal	
PTP-10014E			adapter cable, using with disposable temperature probe; cable material is TPU
PTP-10014F			adapter cable, dual head
PTP-10015A	atom incubator temperatu re probe	<p>skin surface probe</p> 	




PTP-10017A	incubator temperature probe	<p>skin surface probe</p> 		<p>FANEM temperature probe; Reference: 092.059.600</p>
PTP-10018A	<b>Biolight</b>	adult skin surface probe	biolight	
PTP-10018B		adult rectal/esophageal		
PTP-10018C		pediatric skin surface		
PTP-10018D		pediatric rectal/esophageal		
PTP-10018E				adapter cable, using with disposable temperature probe; cable material is TPU




PTP-10019A	creative 25°C @10KΩ		creative	
PTP-10019B		adult skin surface probe		
PTP-10019C		adult rectal/esophageal		
PTP-10019D		pediatric skin surface		
PTP-10019E		pediatric rectal/esophageal		adapter cable, using with disposable temperature probe; cable material is TPU
PTP-10020A	siemens temperature 25°C @6.68KΩ			
PTP-10020B		adult skin surface probe		
PTP-10020C		adult rectal/esophageal		
PTP-10020D		pediatric skin surface		








PTP-10020D		pediatric rectal/esophageal		
PTP-10021A	<b>drager</b> <b>25°C @6.68KΩ</b>	adult skin surface probe		
PTP-10021B		adult rectal/esophageal		
PTP-10021C		pediatric skin surface		
PTP-10021D		pediatric rectal/esophageal		
PTP-10022A	<b>series</b>	adult skin surface probe		
PTP-10022B		adult rectal/esophageal		



PTP-10022C	18B20	pediatric skin surface	
PTP-10022D		pediatric rectal/esophageal	
PTP-10023A	2.25K ohm@25°C GE/ Ohmeda Medical GIRAFFE & GIRAFFE Omni-Bed Incubators and PANDA Warmers	adult skin surface probe	
PTP-10023B		adult rectal/esophageal	
PTP-10023C		pediatric skin surface	
PTP-10023D		pediatric rectal/esophageal	
PTP-10024		disposable temperature probes pads	

PTP-10025	2.25K $\Omega$ @ 25 $^{\circ}\text{C}$	disposable temperature probe;	
PTP-10026A	ningbo david infant warmer; 6 pins	adult skin surface probe	
PTP-10026B		adult rectal/esophageal	
PTP-10026C		pediatric skin surface	
PTP-10026D		pediatric rectal/esophageal	
PTP-10027A	jona	adult skin surface probe	
PTP-10027B		adult rectal/esophageal	

PTP-10027C	arac	pediatric skin surface		
PTP-10027D		pediatric rectal/esophageal		
PTP-10028A	ningbo david infant warmer; 4 pins	adult skin surface probe		
PTP-10028B		adult rectal/esophageal		
PTP-10028C		pediatric skin surface		
PTP-10028D		pediatric rectal/esophageal		
PTP-10029A	re probe,	adult skin surface probe		

PTP-10029B	compatible Zoncure temperaturi new	adult rectal/esophageal			
PTP-10029C		pediatric skin surface			
PTP-10029D		pediatric rectal/esophageal			
PTP-10030A		adult skin surface probe			
PTP-10030B		adult rectal/esophageal			
PTP-10030C		pediatric skin surface			
PTP-10030D		pediatric rectal/esophageal		compatible with Fanem ref: 096.200.600	

PTP-10031A		adult skin surface probe		
PTP-10031B		adult rectal/esophageal		
PTP-10031C		pediatric skin surface		
PTP-10031D		pediatric rectal/esophageal		

